


TEMPERATURE MANAGEMENT BLANKETS


FEATURES & BENEFITS

- **Maximum Mobility:** Moves with the patient from pre-op through intra-op to post-op
- **Cost-Effective:** Only ONE blanket needed for each patient
- **Versatility:** Treats and helps prevent shock, trauma, and hypothermia in a wide variety of settings
- **Quality Materials:** Precision engineered thermal pads, medical grade bi-component non-woven material and proprietary high frequency welding
- **Air-Activated:** No electricity required
- **Warms In Minutes** to 104°F/40°C and maintains temperature for 10+ hours
- **Fully Disposable:** Self-contained system does not circulate particles and one-time use blanket reduces risk of cross contamination
- **US Military Medical Protocol** includes Ready-Heat blankets as part of Tactical Combat Casualty Care since 2006
- **Five Year Shelf Life** and compact storage



Full Body Blanket

Part No.: SB9RH9120
Open Size: 34" x 60" (86 cm x 152 cm)
Description: Blue medical grade non-woven material with 9 oxygen activated engineered thermal pads.



Half Body Blanket

Part No.: SB6RH9120
Open Size: 24" x 34" (60 cm x 86 cm)
Description: Blue medical grade non-woven material with 6 oxygen activated engineered thermal pads.

MANUFACTURED BY TECHTRADE LLC IN THE USA 


Visit our website for more detailed information and IFUs - www.ready-heat.com/IFU
www.ready-heat.com • info@ready-heat.com • 1-866-491-TECH • F: +1 201 706 8130



4 PANEL BLANKET

NSN #: 6532-01-525-4062

Part Number / Per Case: S4RHMD / Case of 8



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 48" / 86 cm x 122 cm


Weight: 1.75 lbs. / 794 gm

Description: Torso Sized Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250 lbs./115 kg). Military Medical Protocol.

6 PANEL BLANKET

NSN #: 6532-01-525-4063

Part Number / Per Case: S6RHLG / Case of 8



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 60" / 86 cm x 152 cm

Weight: 2.5 lbs. / 1.13 kg

Description: Full Body Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250lbs./115 kg). Number 1 choice of US Air Force transports.

1 PANEL BLANKET

NSN #: 6532-01-538-1525

Part Number / Per Case: S1RHSM / Case of 12



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 16" x 18" / 41 cm x 46 cm


Weight: 7 oz. / 198 gm

Description: Mini Blanket. Suitable for injured extremities. Fully adjustable with 2 hook and loop closures. Perfect for First Aid and Survival Kits.

READY-HEAT II BLANKET

NSN #: 6532-01-538-0759

Part Number / Per Case: G12RHZ / Case of 8



Warming Temperature:
100°F / 38°C in 8-10 Minutes

Maintains Heat: Up to 8 Hours


Size: 34" x 48" / 86 cm x 122 cm

Weight: 1.5 lbs. / 794 gm

Description: First Responder Blanket. Used to Treat Shock, Trauma and Pre-Hypothermia. Has a foam core for water resistance and is lightweight and comfortable. Quick to reach operating temperature. Works well at high altitudes.

READY-HEAT II TORSO BLANKET

Part Number / Per Case: G6RH2T / Case of 12



Warming Temperature:
100°F / 38°C in 8-10 Minutes

Maintains Heat: Up to 8 Hours

Size: 34" x 28" / 86 cm x 61 cm

Weight: 14 oz. / 397 gm

Description: First Responder Blanket. Used to treat shock, trauma and pre-hypothermia. Smaller overall footprint both in bag and on patient. All the benefits of the full-sized blanket, G12RHZ.

READY-HEAT PERSONAL USE VEST

Part Number / Per Case: GRHV06A / Case of 12



Warming Temperature:
100°F / 38°C in 10-15 Minutes

Maintains Heat: Up to 8 Hours


Size: Overhead Style with Adjustable Hook and Loop Straps. One Size Fits All.

Weight: 10 oz. / 284 gm

Description: Designed to keep the First Responder warm during long missions. Also great for outdoor activities (Hiking, Skiing, Camping, Hunting, Fishing, Marathons and Mountain Climbing).

INFANT WARMING COCOON

Part Number / Per Case: GIW6C / Case of 12



Warming Temperature:
100°F / 38°C in 15 Minutes

Maintains Heat: Up to 8 Hours


Size: 31.5" x 28" / 80 cm x 71 cm

Weight: 12 oz. / 340 gm

Description: Award Winning Infant Warming Cocoon is designed for the most precious cargo. Six warming pads to cover the vital organs and two adjustable hook and loop fasteners for size adjustment and attached hood for added heat retention.

TEMPERATURE MANAGEMENT FULL BODY

Part Number / Per Case: SB9RH9120 / Case of 10



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 60" / 86 cm x 152 cm

Weight: 18 oz. / 510 gm

Description: Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Prevents hypothermia, increases patient comfort for higher quality of care, fully disposable, self-contained system.

TEMPERATURE MANAGEMENT HALF BODY

Part Number / Per Case: SB6RH9120 / Case of 20



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 24" / 86 cm x 61 cm

Weight: 12 oz. / 350 gm

Description: Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Half-sized blanket can be moved and adjusted for different types of surgeries or access. All the benefits of the full-sized blanket, SB9RH9120.

MANUFACTURED BY TECHTRADE LLC IN THE USA 

Visit our website for more detailed information and IFUs - www.ready-heat.com/IFU
www.ready-heat.com • info@ready-heat.com • 1-866-491-TECH • F: +1 201 706 8130


0459


ISO 13485:2016

