


Monitor Arms

SSM Monitor Arm

Side to Side



Features:

- 180° side to side motion
- Can be mounted directly to any equipment rail or wall channel
- Integrated cable management
- Color options available (white is standard)
- 9" or 12" (23 cm or 30 cm) single and dual arm lengths available
- Swivel arm provides lateral functionality
- Tilt and swivel adjustments
- 30 lbs (13.6 kg) single arm rating
- 15 lbs (6.8 kg) extension arm rating



Tilt up




Level



Tilt down

Tilt and pivot monitor shelf



Fully adjustable tilt and pivot

Monitor Compatibility

PATIENT MONITORS

- | | | |
|-----------------|---------------|-----------|
| • Draeger | • Welch Allyn | • GE |
| • Spacelabs | • Datascope | • Invivo |
| • Fukuda Denshi | • Philips | • Mindray |
| • Nihon Kohden | • Criticare | • Masimo |
| • Datex Ohmeda | • Infinium | |

LCD MONITORS

- | | | |
|-----------|-------------|-------|
| • Toshiba | • Gateway | • NEC |
| • Sony | • Advantech | • LG |
| • Samsung | • DELL | |

Options available for other manufacturers please specify model when ordering

Mounting Options for SSM Models

Channel Mount Options



Channel to **Wall**

- Channel mounts directly to the wall

Part Number: SSM-150



Channel to **Power Column**

- SSM mounts directly to existing monitor channel on power column or flatwall

Part Number: SSM-150



Channel to **Headwall**

- Channel mounts directly to existing rail on headwall
- Channel can be mounted anywhere along length of headwall

Part Number: SSM-150



MRS Channel to **Equipment Rail**

- Channel mounts directly to equipment rails
- Channel can be mounted anywhere along length of rail

Part Number: SSM-150

Direct Rail Mount Options

Benefits of **Direct Rail Mounting**

- Eliminates the need for a monitor channel
- Cleaner, more aesthetically pleasing look
- Increases flexibility for location of monitor
- 40 lbs (18 kg) maximum weight rating (less for some models)
- Compatible with various headwall, flatwall, equipment booms and rail systems



Direct to **Power Column**

- Mounts directly to existing vertical rail, on power column or flatwall


Part Number: SSM-050-V



Direct to **Headwall**

- Mounts directly to existing horizontal rail on headwall

Part Number: SSM-050-A



Direct to **Equipment Rail**

- Mounts directly to equipment rail

Part Number: SSM-LCD-45-ARS

Mounting options are available for both AHM and SSM Arm Models

SSM Models (Side to Side Motion)



Flush Mount

- 4.5" (11.4 cm) profile
- 30 lbs (13.6 kg) single arm rating
- 15 lbs (6.8 kg) dual arm rating
- Back mount devices only
- Quick release auto locking plunger

Part Number: SSM-456



LCD Flush Mount

- 4.5" (11.4 cm) profile
- 30 lbs (13.6 kg) single arm rating
- 15 lbs (6.8 kg) dual arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing

Part Number: SSM-LCD-45

Pivot Arms

 SSM 9" Single (22.9 cm)	 SSM 12" Single (30.5 cm)
13" (33 cm) reach* Part No: SSM-25X	15" (38 cm) reach* Part No: SSM-05X

*30 lbs (13.6 kg) max. per arm

 SSM 9"x9" Double (22.9 cm x 22.9 cm)	 SSM 12"x12" Double (30.5 cm x 30.5 cm)
18" (45.7 cm) reach** Part No: SSM-35X	24" (61 cm) reach** Part No: SSM-15X

**15 lbs (6.8 kg) max. per arm

 LCD SSM 9" Single (22.9 cm)	 LCD SSM 12" Single (30.5 cm)
12" (30.5 cm) reach* Part No: SSM-LCD-25	15" (38 cm) reach* Part No: SSM-LCD-05


*30 lbs (13.6 kg) max. per arm

 LCD SSM 9"x9" Double (22.9 cm x 22.9 cm)	 LCD SSM 12"x12" Double (30.5 cm x 30.5 cm)
19" (48.3 cm) reach** Part No: SSM-LCD-35	24" (61 cm) reach** Part No: SSM-LCD-15

**15 lbs (6.8 kg) max. per arm


SSM Monitor Arm Part Matrix

SSM-XXX-(Z)




LCD SSM Monitor Arm Part Matrix

SSM-LCD-XX-(Z)



Standard Plates Provided:



Monitor Arm Accessories

Monitor Specific Peripheral Mounts Available

Power Bar



- Hospital grade, 4 power outlets
- 120 VAC, 60 Hz

Part Number: PWRBAR4-MON

Power Holder



- Holds most power adapters for PCs, patient monitors and LCDs

Part Number: PWRHLD-MON

Cord Hook



- Manage cables, cords and other hanging wires
- Attach to the bottom of any monitor arm

Part Number (2 hooks): CRDHLDLDR
Part Number (4 hooks): CRDHLDLR2

Bar Code Scanner Holder




- Universal bar code scanner holder

Call for part number.

LCD Handle

- Ergonomic handles can be added to any Amico LCD mount
- Knuckle guard protects walls and user's hands



Part Number: MON-LCD-HANDLE

Cord Wrap



- Ideal for managing wires and cords

Part Number: CORDWRAP

Monitor Channel



- Standard wall channel 34" (48.3 cm)
- Meets seismic requirements
- Custom lengths are welcome

Part Number: MRS-S034

Channel Covers



- Easy way to hide wires
- Available in patterns and colors to match any decor

Call for part number.


Monitor Handle




- Easy grab handle

Part Number: AHM-HAND-00


Wire Baskets



10" W x 10" D x 4" H
(25.4 cm x 25.4 cm x 10.16 cm)
White wire basket
Part No: 101004-MON-BB



10" W x 06" D x 09" H
(25.4 cm x 15.24 cm x 22.86 cm)
White wire basket
Part No: 100609-MON-BB



13" W x 09" D x 07" H
(33.02 cm x 22.86 cm x 17.78 cm)
White wire basket
Part No: 130907-MON-BB



13" W x 09" D x 07" H
(33.02 cm x 22.86 cm x 17.78 cm)
BP cuff basket
Part No: CUFBASK-MON

All baskets can be sterilized

CPU Mounts



- Available Sizes:
- 1.5" - 3" (3.81 cm - 7.62 cm)
Part No: CPUMNT1-MON
- 1.5" - 7" (3.81 cm - 17.78 cm)
Part No: CPUMNTA-MON
- 7" - 9.5" (17.78 cm - 24.13 cm)
Part No: CPUMNT4-MON



- Thin client mount behind LCD
- 10" x 10" x 1.75"
(25.4 cm x 25.4 cm x 4.45 cm)
- 10" x 10" x 3"
(25.4 cm x 25.4 cm x 7.62 cm)



- Enclosed CPU Mount



- Tangent Rugged Mini mount

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more Mounting Solutions.



Distributed By:

